

Protein & Enzyme Lab (BBT 314)

Experiment 2 A:

Determination of the enzyme ALT or SGPT activity in serum by enzymatic method using Bioanalyzer

Background:


Alanine aminotransferase (glutamate pyruvate transaminase) belongs to the group of transaminases, which catalyse the conversion of amino acid to the corresponding α -keto acids via the transfer of amino groups; they also catalyse the reverse process. Although higher activities exist in the liver, minor activity can also be detected in the kidneys, heart, skeletal muscle, pancreas, spleen and lungs. Elevated serum ALT is found in hepatitis, cirrhosis, obstructive jaundice, carcinoma of the liver and chronic alcohol abuse. ALT is only slightly elevated in patients who have an uncomplicated myocardial infarction. Although, both serum aspartate aminotransferase (AST) and ALT found to be elevated whenever disease processes affect liver cell integrity, ALT is more liver specific enzyme. Moreover, elevations of ALT persist longer than elevations of AST activity.

Clinical Significance:

ALT is a cellular enzyme, found in highest concentration in liver and kidney. High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism, its better application is in diagnosis of the disease of the liver.

Principle:

Alanine aminotransferase (glutamate pyruvate transaminase) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and pyruvate. The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH.



The rate of decrease in concentration of NADH measured photometrically at 340nm is proportional to the catalytic concentration of ALT present in the sample.

Reagents:

R1 (buffer)	Tris pH 7.8	100mmol/l
	L-Alanine	500mmol/l
	Lactate dehydrogenase (LDH)	1200U/l
R2(substrate)	NADH	0.18mmol/l
	α -ketoglutarate	15mmol/l

Specimen: Serum

Specimen Preparation:

Separate serum from cells by centrifugation. ALT activity is stable at 2-8°C for 7 days.

Reagent preparation:

Working reagent (WR): mix 1 volume of R2 with 4 volume of R1. Reagent is stable for 21 days at 2-8°C and 72 hours at room temperature.

Working reagent (WR)	
R1	4 volume
R2	1 volume

Procedure:

1. Set the spectrophotometer at 340 nm and adjust to zero with dH₂O. This is the blank reading.
2. Prepare a tube as follows,

WR	1.0 ml
Sample	100 µl

Mix and incubate for 1 minute.

3. Transfer about 200µl to a cuvette and read the absorbance of the sample. This is the absorbance reading at 0 minute.
4. Start the stopwatch and read the absorbance at 1 minute interval for 3 minutes at 340 nm.

Reading no.	Time (minute)	Absorbance
1	1	
2	2	
3	3	
4	4	

Calculation:

Calculate the difference of absorbance and the average absorbance difference per minute ($\Delta A/\text{min}$).

Difference of absorbance (ΔA)		
$\Delta A1$	absorbance at 1 min-absorbance at 2 min	
$\Delta A2$	absorbance at 2 min-absorbance at 3 min	
$\Delta A3$	absorbance at 3 min-absorbance at 4 min	

Average Absorbance difference per minute ($\Delta A/\text{min}$) = $(\Delta A1 + \Delta A2 + \Delta A3)/3$

Write down the calculation of ALAT/ALT/GPT concentration by using given formula.

$$\text{GPT/ALT} = \Delta A/\text{min} \times 1750 \text{ U/L}$$

Reference value (At 37°C):

Men: up to 40U/L **Women:** up to 32U/L

Experiment 2 B

Determination of the enzyme AST or GOT activity in serum by enzymatic method using Bioanalyzer

Background:


Aspartate aminotransferase (glutamate oxaloacetate transaminase) belongs to the group of transaminases, which catalyze the conversion of amino acids to the corresponding α -keto acids via the transfer of amino groups; they also catalyze the reverse process. AST is commonly found in human tissue. Although heart muscle is found to have the most activity of the enzyme, significant activity has also been seen in the brain, liver, gastric mucosa, adipose tissue, skeletal muscle, and kidneys. AST is present in both cytoplasm and mitochondria of cells. In cases with mild tissue injury, the predominant form of AST is that from the cytoplasm, with a smaller amount coming from the mitochondria. Severe tissue damage results in more of the mitochondrial enzyme being released. Elevated levels of transaminases are indicative of myocardial infarction, hepatopathies, muscular dystrophy and damage to the internal organs. Increased levels of AST however are generally a result of viral or toxic hepatitis and obstructive jaundice. Following a myocardial infarction, serum levels of AST are elevated and reach a peak 48-60 hours after onset.

Clinical Significance:

The AST is a cellular enzyme, is found in highest concentration in heart muscle, the cells of the liver, the cells of the skeletal muscle and in smaller amounts in other tissues. Although an elevated level of AST in the serum is not specific of the hepatic disease, is used mainly to diagnostic and to verify the course of this disease with other enzymes like ALT and ALP. Also it is used to control the patients after myocardial infarction, in skeletal muscle disease and other^{1, 4, 5}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

Principle:

Aspartate aminotransferase (AST) formerly called glutamate oxaloacetate (GOT) catalyses the reversible transfer of an amino group from aspartate to alpha-ketoglutarate forming glutamate and oxalacetate. The oxalacetate produced is reduced to malate by malate dehydrogenase (MDH) and NADH:



Reagents:

R1 (buffer)	Tris pH 7.8	100mmol/l
	L-Alanine	500mmol/l
	Lactate dehydrogenase (LDH)	1200U/l
R2 (substrate)	NADH	0.18mmol/l
	α -ketoglutarate	15mmol/l

Specimen: Serum

Specimen Preparation:

Separate serum from cells by centrifugation. AST activity is stable at 2-8°C for 7 days.

Reagent preparation:

Working reagent (WR): mix 1 volume of R2 with 4 volume of R1. Reagent is stable for 21 days at 2-8°C and 72 hours at room temperature (15-25°C).

Working reagent (WR)	
R1	4 volume
R2	1 volume

Procedure:

1. Set the spectrophotometer at 340 nm and adjust to zero with dH₂O. This is the blank reading.
2. Prepare a tube as follows,

WR	1.0 ml
Sample	100 μ l

3. Mix and incubate for 1 minute.
4. Transfer about 200 μ l to a cuvette and read the absorbance of the sample. This is the absorbance reading at 0 minute.
5. Start the stopwatch and read the absorbance at 1 minute interval for 3 minutes at 340 nm.

Reading no.	Time (minute)	Absorbance
1	1	
2	2	
3	3	
4	4	

Calculation:

Calculate the difference of absorbance and the average absorbance difference per minute (ΔA /min).

Difference of absorbance (ΔA)		
ΔA_1	absorbance at 1 min-absorbance at 2 min	
ΔA_2	absorbance at 2 min-absorbance at 3 min	
ΔA_3	absorbance at 3 min-absorbance at 4 min	

$$\text{Average Absorbance difference per minute } (\Delta A/\text{min}) = (\Delta A_1 + \Delta A_2 + \Delta A_3)/3$$

CALCULATIONS

$$\Delta A/\text{min} \times 1750 = \text{U/L of AST}$$

Reference value (At 37°C):

Men: up to 38 U/L **Women:** up to 31 U/L